



Standardized result verification for faster, consistent turn-around time

With a shrinking pool of qualified personnel to produce results, you need a system that will eliminate mundane tasks and allow your scarce resources, medical technologists, to focus on the 'true' exceptions. Whether you are a reference or hospital laboratory, have one LIS or five, autoverification from Data Innovations' middleware delivers the most flexible rules-based decision processing to verify your results with maximum consistency and accuracy.

Notifier Alert

With Data Innovations' unique, powerful Notifier application, you tie your event notifications – and the recipients – directly to the rules you choose. (For example, getting an email or text message the moment an elevated Troponin result is produced.)

Is your lab **IM**powered™ ?



Scrutinizes results with maximum consistency and accuracy

From basic to complex, the rules you write are there to help you control your workflow and improve quality results and services you want to provide. That's why we deliver total flexibility. You may create the rules needed to implement autoverification to cover all lab disciplines in all locations, while maintaining management of the rules from within the lab.

- Rules-based decision processing via graphical interface with drop down menus
- Easy to develop rules to support serum indices, real-time delta checking, QC integration, eGFR, reflex and anti-reflex testing and more
- Offline rules test environment lets you simulate information from the instrument and/or LIS to verify your rules logic
- Print rules testing for regulatory compliance
- User-defined, customized results review screens support graphics, previous results, filtering and color coding using rules logic