



Access to critical information for improving analysis and operations

No matter how many sites, instruments or customers you need to serve, Data Innovations' middleware conveniently consolidates all your analytical data in one location for easy retrieval for analysis and reporting.

Notifier Alert

With Data Innovations' unique, powerful Notifier application, you tie your event notifications – and the recipients – directly to the rules you choose. (For example, if you are looking for a sample that meets certain criteria for a study.)

Is your lab **IM**powered™ ?



Built-in flexibility to track, report, and review your data

You can collect data by push (data collection module) or pull (ODBC access) methods. Query the database and analyze data using virtually any third-party application. We've built in all the flexibility you need.

- Ad hoc report generation with any ODBC compliant application such as Microsoft® Excel or Crystal Reports®
- Automated generation of information and delivery via email or file server on a daily, weekly, monthly basis
- Logical names for database record fields for easier data mining
- Eases mundane tasks such as correlations for instrument or lot comparisons, determining moving averages, TAT calculations, and more