



Automate your hematology workarea in days, not months

The motivation for autoverification is evident; increased financial constraints, increasing volumes of testing, increased regulatory constraints, demands for decreased response time, and higher expectations on quality.

The motivation is there, the spare time to do it is not. IM Solutions - Hematology from Data Innovations provides a comprehensive, turnkey middleware solution for your hematology workarea to optimize and automate result verification in days, not months.

Is your lab **IM**powered™ ?



Scrutinize results with maximum consistency and accuracy

A comprehensive hematology rule set for autoverification eliminates many mundane tasks and reduces staff fatigue. The “true” exceptions are presented on one concise screen that consolidates information from the LIS, the instrument(s), the microscope, and previous patient information to provide a holistic approach for results review.

- Turnkey middleware solution with hardware, software, onsite installation, onsite training, and support included.
- Support for multiple hematology instruments, LIS connectivity and multiple, concurrent users.
- A complete hematology rule set for autoverification that includes the International Society for Laboratory Hematology (ISLH) Consensus Guidelines.
- Pre-defined, customizable results review screens that support instrument graphics, previous results, comments, and integrated cell counter to perform cell morphology and manual differential counting.
- Different configurations are optimized to work with the leading hematology instruments.