



Taking Quality Control to the next level

The primary service of your laboratory is providing actionable healthcare information. With integrated QC, you will strengthen your confidence in your laboratory's precision and reliability of its results.

Notifier Alert

With Data Innovations' unique, powerful Notifier application, pop-up notifications deliver critical flagged information to operators immediately so they can take appropriate action to maintain your quality standards. (For example, when a 1-3s rule violation occurs or QC has not been run in 8 hours for a particular instrument.)

Is your lab **IM**powered™ ?



Unique bi-directional functionality strengthens your QC program

We are the only middleware company to offer a bi-directional interface with Bio-Rad's QC OnCall™ and Unity Real Time™. That's because we recognize the vital importance of integrating QC into results production. Nowhere is this more critical to your business than in quality control procedures to protect your physicians and their patients.

- Automatically and immediately send QC results to QConCall™ / Unity Real Time™
- Integrate QC rules evaluation within Instrument Manager autoverification rules allowing you to alter results production immediately
- On-demand interlaboratory comparison reports for real-time troubleshooting
- Web-based, peer grouping via Bio-Rad's QC Net