

**Thermo Scientific  
Data Innovations and LabLink xL  
Connectivity Solutions**



**With the option of direct connectivity to Data Innovations® Instrument Manager™, laboratories can now benefit from fully-automated transfer of QC data into the LabLink xL QAP system.**

**QC results are uploaded within seconds to LabLink xL Web for immediate review on any PC with access to the Internet. Peer group comparison and report production is available instantaneously through the 100% live LabLink xL database. There are no monthly cut-off dates for data submission - your QC data is never late.**



## Key System Features:

- Fully automatic transfer of QC data from laboratory instrumentation directly into LabLink xL Web
- Instantaneous access to worldwide peer group comparison
- Requires simple installation of Instrument Manager LabLink xL driver
- Access to QC results from any PC with access to the Internet
- Password protected security - data is held securely

## LabLink xL

Used in conjunction with MAS controls, LabLink xL fulfills various regulatory requirements for QC tracking, including peer data review.

- Get an external perspective on your internal day-to-day QC
- Troubleshoot variances immediately with visibility to a world-wide peer group of users
- Assess quickly and easily whether shifts in QC recoveries are unique to your laboratory
- Reduce troubleshooting time and quickly release patient results

## LabLink Extra

LabLink Extra works in unison with the LabLink xL Quality Assurance Program to provide a gateway for up to the minute support data for the QC products manufactured by Thermo Fisher Scientific and the LabLink xL OAP system.

- Access to online QC package insert documents and updates
- Quality Control handling and enhanced performance recommendations
- LabLink xL QC data file generation guidelines (LIS and Instrument Options)
- Country specific technical support contact information available on demand

## Ordering Information

For more information on MAS Quality Controls or LabLink xL and Extra, please contact the Network Sales group at 1-800-232-3342, option 4 or email to: [network.sales@thermofisher.com](mailto:network.sales@thermofisher.com)

### Diagnostics

**U.S.A.**  
46360 Fremont Blvd.  
Fremont, CA 94538 USA  
Tel. 800-232-3342  
Fax: 510-979-5002

**Canada**  
Tel. 905-286-4290  
[contact@diagnostix.ca](mailto:contact@diagnostix.ca)

**Germany**  
Tel. + 49 851 88689 0  
[microgenics.DE.info@thermofisher.com](mailto:microgenics.DE.info@thermofisher.com)

**Australia**  
Tel. + 61 2 9649 9599  
[info@microgenics.com.au](mailto:info@microgenics.com.au)

**Spain**  
Tel. + 34 93589 8338  
[microgenics.ES.info@thermofisher.com](mailto:microgenics.ES.info@thermofisher.com)

**United Kingdom**  
Tel. + 44 1727 853151  
[microgenics.UK.info@thermofisher.com](mailto:microgenics.UK.info@thermofisher.com)